


ROLLEI RXN FIX NEUTRAL

Rollei RXN Fix Neutral is a high-performance, low-odor fixing bath with neutral pH value. RXN is specially designed for the archival processing of baryta papers but is also universally suitable for processing RC papers and black and white films.



Dilution 1 + 4 :



800ml water + 200ml RXN = 1l working solution
Temp. 20°C / 68°F

VERDÜNNUNG	BARYTPAPIER	PE / RC PAPIER	FILM
1+4	2 min.	1 min.	6 min. *
1+7	2-3 min.	1-2 min.	-
1+9	-	2-3 min.	-

* Stop bath required, except for pyro developers (e.g. PMK, Tanol, etc.)

Hint: RXN is reducing the general watering time of FB papers. Please pre-wash shortly paper from the developer in extra clear water before putting in stopping bath. This will reduce lab odor disturbances additionally.

AVAILABLE SIZES



Due to the chemical formula no safety advice is necessary – nevertheless we recommend a responsible handling of the concentrate. Wear suitable protective gloves and clothing. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. **Keep out of the reach of children.**

GENERAL NOTES:

CHEMISTRY | DILUTION | TIME | INTERVAL:
